



# Rapid quantitative analysis of daidzein and genistein in soybeans (*Glycine max* (L). Merr.) using FTIR spectroscopy and multivariate calibration

Etty Sulistyowati<sup>1,2</sup>, Sudibyo Martono<sup>2</sup>, Sugeng Riyanto<sup>2</sup>, Endang Lukitaningsih<sup>2</sup>, Abdul Rohman<sup>2,3\*</sup>

<sup>1</sup>Pharmaceutical Studies, College of Pharmaceutical Science, Yayasan Pharmasi Semarang, Semarang 50192, Indonesia.

<sup>2</sup>Department of Pharmaceutical Chemistry, Faculty of Pharmacy, Universitas Gadjah Mada, Yogyakarta, 55281, Indonesia.

<sup>3</sup>Institute for Halal Industry & System (IHIS), Universitas Gadjah Mada, Yogyakarta, 55281, Indonesia.

## ARTICLE INFO

Received on: 08/05/2020  
Accepted on: 25/07/2020  
Available online: 05/11/2020

### Key words:

Daidzein, FTIR, genistein, PLS, PCR, soybean.

## ABSTRACT

Soybean with the scientific name of *Glycine max* (L.) Merr. is nutritious vegetable food sources in Indonesia, and it may be processed into various kinds of food products. Soybeans are a rich source of isoflavones, rather than nuts and meat products. The research was aimed at quantifying daidzein (DN) and genistein (GN) in various soybean varieties applying Fourier transform Infrared (FTIR) spectroscopy coupled with two multivariate calibrations [principal component regression (PCR) and partial least square (PLS)]. Before being used as variables during quantification processes, some pretreatments of FTIR spectra including the selection of wavenumbers and spectral derivative (high-order derivatives) were carried out to get the best prediction models for the correlation between the actual values of DN and GN. The actual values of DN and GN were determined with High-performance liquid chromatography (HPLC). The results revealed that PLS provides a better modeling than PCR for the relationship between actual values and FTIR spectroscopy-predicted values. The first derivative spectra at wavenumbers of 3,600–2,800 and 1,500–780 using PLS regression were preferred for the quantification of DN and GN in soybean. PLS calibration model yielded  $R^2$  for the relationship between HPLC-based actual values ( $x$ -axis) and FTIR-predicted values ( $y$ -axis) of DN and GN, which were 0.9947 and 0.9900, respectively. FTIR spectroscopy method-combined PLS regression provides fast and acceptable results as indicated by high accuracy and precision during quantification of DN and GN in soybean.

## INTRODUCTION

Soybean [*Glycine max* (L.) Merr.] is a popular health nutritious vegetable food widely consumed in Indonesian. The economic value of soybean is very important for Indonesian society because various processed soy-based food products become their daily consumption. Isoflavone (aglycone form of glycosides) such as daidzein (DN) and genistein (GN) is reported to be more responsible for biological activities than glycoside form because aglycones had a good biological activity and are

also easy to be absorbed (Izumi *et al.*, 2000). DN and GN have potential protective effects on mammary cancers (Dhananjaya *et al.*, 2012), estrogenic effect (Islam *et al.*, 2008), and protective action against osteoporosis (Morabito *et al.*, 2002).

High-performance liquid chromatography (HPLC) is an analytical technique widely used for the analysis of DN and GN in soybean samples or soybean-based food products (Hong *et al.*, 2011; Magiera and Sobik, 2017; Shao *et al.*, 2011; Yatsu *et al.*, 2016). However, an HPLC analysis needs complex sample preparation and skillful analysis; therefore, fast and reliable analytical techniques based on Fourier transform Infrared (FTIR) spectroscopy methods combined with chemometrics are developed for the analysis of DN and GN in soybeans (Mulsow *et al.*, 2015). FTIR spectroscopy is based on the interaction between infrared radiation and samples to get specific peaks corresponding to the absorption of functional groups present in the analyzed samples (Rohman, 2012).

### \*Corresponding Author

Abdul Rohman, Department of Pharmaceutical Chemistry, Faculty of Pharmacy, Universitas Gadjah Mada, Yogyakarta, 55281, Indonesia.  
E-mail: [abdulrohmanugm@gmail.com](mailto:abdulrohmanugm@gmail.com)

FTIR spectroscopy coupled with multivariate data analysis (chemometrics) of multivariate calibrations can be used for the qualitative and quantitative analysis of compounds of natural ingredients (Rohman *et al.*, 2014). FTIR spectroscopy-based techniques have been used for the quantification of curcuminoids in *Curcuma* syrups (Prabaningdyah *et al.*, 2018) and analysis of curcumin and demethoxycurcumin in some extracts of turmeric (*Curcuma longa* L.) and java turmeric (*Curcuma xanthorrhiza*) (Lestari *et al.*, 2017; Rohman *et al.*, 2015). Based on this application, the correlation between FTIR spectroscopic and HPLC methods emerged very flexible for the analysis of DN and GN. Using a literature study, there are no reports regarding the employment of FTIR spectroscopy-based techniques for the analysis of DN and GN. For this reason, this study was aimed to optimize and develop FTIR spectroscopy coupled with partial least square (PLS) and principal component regression (PCR) for the analysis of DN and GN. As actual values of DN and GN in soybeans, the validated HPLC method was used.

## MATERIALS AND METHODS

The reference standards of DN and GN were purchased from Sigma (Aldrich, St. Louis, MO). HPLC instrument was used for the quantitative analysis of DN and GN. Soybean samples were obtained from Balitkabi (Balai Penelitian Tanaman Aneka Kacang dan Umbi), Malang, East Java, Indonesia.

### Preparation of soybean samples

Approximately 2.5 g of soybean samples were extracted with 25 ml of water–methanol (50:50 volume/volume) at room temperature for 3 days with immediate shaking. The mixture was then filtered, and then, 1.0 ml of sample solution was diluted with 9.0 ml of water–methanol (50:50 volume/volume). This solution was subjected to filtration using 0.45- $\mu\text{m}$  polyvinylidene fluoride or polyvinylidene difluoride and finally introduced into HPLC system.

### HPLC analytical conditions

Reserved-phase HPLC for the quantification of DN and GN in soybean samples was performed on HPLC waters, using column Sun Fire <sup>TM</sup>C-18 (150  $\times$  4.6 mm, with internal diameter 5  $\mu\text{m}$ ), equipped with guard column Waters Symmetry <sup>TM</sup>C18 (20  $\times$  4.4 mm, with an internal diameter of 5  $\mu\text{m}$ ). The composition of mobile phase consisted of a mixture of methanol and 0.1% acetic acid (53:47 volume/volume) delivered isocratically using a flow rate of 1.0 ml/minute. The volume of injection was 10  $\mu\text{l}$ . For the preparation of standard solutions, 10 mg of reference standards were accurately weighed and then diluted with mobile phase to get stock solution of DN and GN (100  $\mu\text{g}/\text{ml}$ ). Aliquots of the solution (2, 3, 4, 5, 6, 7, 8, 9, and 10  $\mu\text{g}/\text{ml}$ ) were prepared and injected into the HPLC equipment. The solutes containing DN and GN were detected by using photodiode array at 254 nm. HPLC running time was 20 minutes with retention times ( $t_r$ ) of DN and GN which were 6.342 and 10.088, respectively.

### FTIR spectroscopy analysis

The analysis of samples using FTIR spectroscopy was performed according to Wulandari *et al.* (2018). The powdered samples were directly placed on an attenuated horizontal total

reflectance (Smart iTR<sup>TM</sup>) accessory. The specification of measurement is as follows:

Software:	OMNIC ver. 9.7
Wavenumbers of scanning:	4,000–650 $\text{cm}^{-1}$
Number of scanning:	32 scans
Spectral resolution:	8 $\text{cm}^{-1}$
Spectral background:	air (environment) spectrum

### Model prediction of DN and GN

The prediction of DN and GN was facilitated by multivariate calibrations of partial least-square regression (PLSR) and PCR. The levels of DN and GN obtained during HPLC were used as actual values, and then, the levels of DN and GN were predicted using the variables of absorbance values at the region of optimized wavenumbers. The correlation models for actual values and predicted values were regressed using PLSR.

### Data analysis

The software of TQ Analyst from Thermo Fisher Scientific Inc. (Madison, WI) was used for data analysis applying multivariate calibrations including PLS and PCR. Statistical parameters used are as follows:

Coefficient of correlation (R) for correlation between HPLC actual values and FTIR predicted values which indicated the accuracy of the model.


The errors of root mean square error of calibration (RMSEC) and error of prediction (RMSEP) indicating the precision of the model (Miller and Miller, 2005).

## RESULTS AND DISCUSSION


Using two models of multivariate calibrations, PLS and PCR were employed for the correlation between HPLC actual values of DN and GN as determined and FTIR predicted values using variables of absorbance values at certain FTIR spectra regions. The HPLC chromatogram obtained during the separation of DN and GN is shown in Figure 1. HPLC conditions used can provide a good system because it has an asymmetry value of 1.2. The relative standard deviation (RSD) values for peak area and peak height were obtained to meet the requirements of maximum RSD value, namely,  $\leq 2\%$ . In addition, HPLC used was validated previously by Sulistiyowati *et al.* (2019).

FTIR spectra at the wavenumbers of 4,000–650  $\text{cm}^{-1}$  of soybean pulverized containing DN were overlaid as shown in Figure 2. Each peak at specific wavenumbers ( $1/\lambda$ ) in these spectra could be attributed by functional groups present in soybean samples. A wide peak at  $1/\lambda$  3,272  $\text{cm}^{-1}$  was the vibration of -OH stretching, associated with hydrogen bond of -OH. Two peaks at  $1/\lambda$  2,925 and 2,855  $\text{cm}^{-1}$  were coming from vibrations of  $\text{CH}_3$ - and  $\text{CH}_2$ - functional groups in stretching modes. A peak at  $1/\lambda$  1,743  $\text{cm}^{-1}$  originated from carbonyl (C=O) groups (Rohman *et al.*, 2015). The presence of these functional groups indicated the presence of DN and GN in soybean samples.

For quantitative analysis using FTIR spectroscopy, some optimizations were performed by selecting the FTIR spectral treatment and derivatization treatments. Based on the optimization applying wavenumbers' region and FTIR spectral treatments (normal, first, and second derivatives) and relying



**Figure 1.** Chromatogram by HPLC with the SunFire™C-18 (150 × 4.6 mm, 5 μm): (a) chromatogram of soybeans extract containing DN and GN and (b) chromatogram spiking performed using 1 ppm of DN and GN standard.



**Figure 2.** The overlay of FTIR spectra of soybean pulverized containing DN and scanned at midinfrared region (4,000–650 cm<sup>-1</sup>).

on the highest coefficient of correlation ( $R$ -value) and lowest values of RMSEC and RMSEP, DN and GN were subjected to quantification using the first derivative spectra at combined  $1/\lambda$  of 3,600–2,800 cm<sup>-1</sup> and 1,500–780 cm<sup>-1</sup> with 10 factors. The selection of the wavenumbers of samples is important because the measurements with not suitable wavenumber regions or less informative are able to reduce the performance of PLS/PCR modeling (El-Gindy *et al.*, 2006). Tables 1 and 2 show the performance of multivariate calibrations of PLS and PCR models applying either normal or derivative spectra (first and second) for the quantitative analysis of DN and GN in soybeans together

with some statistical performances, consisting of the number of factors,  $R$ -values, and error values (RMSEC and RMSEP).

The values of  $R$  obtained for the determination of DN using PLS were 0.9943 (calibration model) and 0.9961 (validation model), whereas the error values were 0.865% (calibration) and 0.994% (validation), respectively (Fig. 3). The  $R$ -values for the determination of GN were 0.9882 (calibration) and 0.9936 (validation) with the error values of 0.874% (calibration) and 1.01% (validation), respectively (Fig. 4). Figure 5 shows the correlation between HPLC actual results ( $x$ -axis) and FTIR calculated values ( $y$ -axis) in validation samples which resulted

**Table 1.** The characteristic performances of PLS for modeling the correlation between HPLC actual values and FTIR predicted values of DN and GN.

Wavenumber (cm <sup>-1</sup> )	Spectral treatment	Number of factors	DN				GN			
			Calibration		Validation		Calibration		Validation	
			RMSEC	R	RMSEP	R	RMSEC	R	RMSEP	R
4,000–650	Normal	10	2.74	0.9416	3.98	0.9349	4.13	0.6911	5.25	0.7318
	Derivative 1	6	3.55	0.8997	3.62	0.9369	2.83	0.8691	4.15	0.8347
	Derivative 2	7	2.67	0.9447	6.17	0.7962	2.07	0.9324	4.91	0.7673
3,600–650	Normal	10	2.67	0.9446	3.54	0.9508	4.1	0.6973	5.31	0.7309
	Derivative 1	6	3.35	0.9115	3.23	0.9524	2.83	0.8684	4.22	0.8283
	Derivative 2	7	2.46	0.9532	6.02	0.8043	0.607	0.9943	4	0.8541
3,600–780	Normal	9	2.83	0.9377	3.77	0.9364	1.84	0.9466	3.22	0.9231
	Derivative 1	10	1.22	0.9887	1.44	0.9905	2.88	0.8636	3.75	0.8739
	Derivative 2	10	0.881	0.9941	4.22	0.9073	0.579	0.9949	3.38	0.8967
1,789–780	Normal	9	2.87	0.9358	3.61	0.9491	2.07	0.9322	2.7	0.9477
	Derivative 1	10	1.31	0.9869	1.68	0.9855	1.6	0.96	1.9	0.9762
	Derivative 2	9	1.66	0.9789	3.54	0.9334	3.14	0.8353	3.76	0.8899
1,500–780	Normal	10	2.52	0.9507	3.68	0.9462	2.32	0.9137	3.13	0.9159
	Derivative 1	9	1.12	0.9905	1.43	0.9902	0.875	0.9882	1.12	0.9923
	Derivative 2	8	1.66	0.9789	3.28	0.9493	1.33	0.9724	2.36	0.9611
3,600–2,800 and 1,789–780	Normal	9	3.2	0.9196	3.18	0.9543	4.09	0.6991	5.23	0.7647
	Derivative 1	8	0.86	0.9734	1.57	0.9873	2.45	0.9031	2.72	0.9522
	Derivative 2	7	2.48	0.9524	3.99	0.9163	2.78	0.8737	4.23	0.832
3,600–2,800 and 1,500–780	Normal	9	3.09	0.9254	3.27	0.9507	2.29	0.9166	3.01	0.9226
	Derivative 1	10	0.865	0.9943	0.994	0.9961	2.27	0.9179	2.53	0.9562
	Derivative 2	5	3.78	0.8855	4.96	0.8662	2.33	0.9135	3.87	0.861
3,000–2,800 and 1,789–780	Normal	6	4.52	0.8314	4.96	0.8863	2.43	0.9052	3.14	0.9207
	Derivative 1	10	1.26	0.988	1.25	0.9919	1.02	0.9839	1.57	0.9794
	Derivative 2	7	2.58	0.9486	3.63	0.9315	3.16	0.833	3.77	0.8939
3,000–2,800 and 1,500–780	Normal	7	3.77	0.8867	4.28	0.9191	2.49	0.9002	3.19	0.9155
	Derivative 1	10	0.996	0.9925	1.26	0.9921	0.874	0.9882	1.01	0.9936
	Derivative 2	4	4.22	0.855	4.84	0.8705	2.76	0.8756	3.7	0.8726

R = coefficient of correlation; RMSEC = root mean square error of calibration; RMSEP = root mean square error of prediction.



**Table 2.** The performance of multivariate calibrations of PCR for modeling HPLC actual values and FTIR predicted values of DN and GN.

Wavenumber (cm <sup>-1</sup> )	Spectral treatment	Number of factors	DN				GN			
			Calibration		Validation		Calibration		Validation	
			RMSEC	R	RMSEP	R	RMSEC	R	RMSEP	R
4,000–650	Normal	10	6.18	0.6504	7.32	0.7253	4.28	0.6634	5.42	0.7366
	Derivative 1	10	5.63	0.7225	6.2	0.7827	3.78	0.7503	5.19	0.7221
	Derivative 2	10	6.98	0.514	8.47	0.5225	4.62	0.5891	6.31	0.5422
3,600–650	Normal	10	6.2	0.6484	7.38	0.7183	4.29	0.6603	5.52	0.7219
	Derivative 1	10	5.52	0.7351	5.85	0.81	3.73	0.7582	4.86	0.7608
	Derivative 2	10	7.06	0.4991	8.48	0.514	4.68	0.5749	6.64	0.469
3,600–780	Normal	10	6.17	0.653	7.22	0.7402	4.27	0.6651	5.41	0.7427
	Derivative 1	10	5.36	0.7527	5.37	0.8483	3.77	0.7523	4.22	0.8347
	Derivative 2	10	6.7	0.5687	7.59	0.6503	4.4	0.638	5.32	0.7226
1,789–780	Normal	10	5.64	0.7214	6.6	0.7567	3.6	0.7771	4.5	0.8226
	Derivative 1	10	5.34	0.7544	5.84	0.8266	3.63	0.7721	3.94	0.8662
	Derivative 2	10	6.09	0.6642	7.31	0.6792	3.82	0.7437	4.54	0.8113
1,500–780	Normal	10	4.78	0.8092	6.09	0.8111	2.81	0.871	3.67	0.8805
	Derivative 1	10	5.23	0.7662	6.21	0.7858	3.7	0.7615	4.69	0.7839
	Derivative 2	10	5.25	0.7647	5.54	0.8306	3.72	0.7596	4.65	0.7976


Continued

Wavenumber (cm <sup>-1</sup> )	Spectral treatment	Number of factors	DN				GN			
			Calibration		Validation		Calibration		Validation	
			RMSEC	R	RMSEP	R	RMSEC	R	RMSEP	R
3,600–2,800 and 1,789–780	Normal	10	6.19	0.6494	7.44	0.7381	4.28	0.6625	5.54	0.7238
	Derivative 1	10	5.33	0.7559	5.82	0.8268	3.57	0.7812	4.1	0.8531
	Derivative 2	10	6.34	0.6274	6.75	0.7501	4.12	0.6937	4.88	0.7786
3,600–2,800 and 1,500–780	Normal	10	5.67	0.7172	6.18	0.8943	4.15	0.687	5.02	0.8455
	Derivative 1	10	5.47	0.7404	6.47	0.7658	3.77	0.7514	4.94	0.7656
	Derivative 2	10	5.58	0.7281	5.29	0.871	3.97	0.7189	5.05	0.7548
3,000–2,800 and 1,789–780	Normal	10	6.07	0.6664	7.39	0.7081	4.29	0.6616	5.84	0.661
	Derivative 1	10	5.32	0.7573	5.75	0.8339	3.61	0.7747	3.92	0.8687
	Derivative 2	10	6.29	0.6357	6.6	0.7616	4.11	0.6952	4.76	0.7918
3,000–2,800 and 1,500–780	Normal	10	4.95	0.7945	6.49	0.7982	3.11	0.8386	4.38	0.8355
	Derivative 1	10	5.06	0.7837	5.38	0.8546	3.56	0.7828	4.24	0.8355
	Derivative 2	10	5.68	0.7166	5.59	0.8425	4.06	0.7047	4.98	0.7621


R = coefficient of correlation; RMSEC = root mean square error of calibration; RMSEP = root mean square error of prediction.



**Figure 3.** PLS model for the relationship between HPLC actual values (x-axis) and FTIR calculated values of DN. (a) Calibration model; (b) residual analysis.



**Figure 4.** PLS model for the relationship between HPLC actual values (x-axis) and FTIR calculated values of GN. (a) Calibration model; (b) residual analysis.



**Figure 5.** Scatter plot for the relationship between HPLC actual values of DN and GN and FTIR predicted values. (a) Daidzein; (b) GN.  $R^2$  = coefficient of determination.

in  $R^2 > 0.99$ . This indicated that PLS with appropriate spectral treatments is accurate and precise methods for the quantification of DN and GN in soybeans (Sim *et al.*, 2004).

## CONCLUSION

HPLC was used successfully for the analysis of DN and GN in soybean samples, and the concentrations obtained were used for actual values during the prediction of DN and GN in soybean samples. FTIR spectroscopy using the first derivative spectra at combined  $1/\lambda$  of 3,600–2,800 and 1,500–780  $\text{cm}^{-1}$  coupled PLS model could be used as a rapid and reliable method for the quantitation of DN and GN in soybeans. FTIR spectroscopy offered the reliable method for the analysis of DN and GN in soybeans with the main advantages of being rapid, minimum sample preparation, and being environmental friendly. The developed method could be prolonged to be used in the analysis of soybean-based food products such as soybean milk and soy sauce.

## ACKNOWLEDGMENTS

The authors acknowledged the Ministry of Research and Higher Education, Republic of Indonesia, for financial support during this study.

## CONFLICT OF INTEREST

Authors declared that they do not have any conflicts of interest.



**FUNDING**

None.

**REFERENCES**

- Dhananjaya K, Sibi G, Mallesha H, Ravikumar KR, Awasthi S. In silico studies of daidzein and genistein with human estrogen receptor  $\alpha$ . *Asian Pac J Trop Biomed*, 2012; 2:1747–53.
- El-Gindy A, Emara S, Mostafa A. Application and validation of chemometrics-assisted spectrophotometry and liquid chromatography for the simultaneous determination of six-component pharmaceuticals. *J Pharm Biomed Anal*, 2006; 41(2):421–30.
- Hong JL, Qin XY, Shu P, Wang Q, Zhou ZF, Wang GK. Comparative study of isoflavones in wild and cultivated soybeans as well as bean products by high-performance liquid chromatography coupled with mass spectrometry and chemometric techniques. *Eur Food Res Technol*, 2011; 233:869–80.
- Islam, F., Sparkes, C., Roodenrys, S., dan Astheimer, L., 2008. Short-term changes in endogenous estrogen levels and consumption of soy isoflavones affect working and verbal memory in young adult females. *Nutr Neurosci*, 11(6): 251–62.
- Izumi T, Piskula MK, Osawa S, Obata A, Tobe K, Saito M. Soy isoflavone aglycones are absorbed faster and in higher amounts than their glucosides in humans. *J Nutrition*, 2000; 130:1695–9.
- Lestari HP, Martono S, Wulandari R, Rohman A. Simultaneous analysis of curcumin and demethoxycurcumin in *Curcuma xanthorrhiza* using FTIR spectroscopy and chemometrics. *Int Food Res J*, 2017; 24:2097–101.
- Magiera S, Sobik A. Ionic liquid-based ultrasound-assisted extraction coupled with liquid chromatography to determine isoflavones in soy foods. *J Food Compos Anal*, 2007; 57:94–101.
- Miller JN, Miller JC. *Statistics and chemometrics for analytical chemistry*. 5th edition, Pearson Education Limited, Edinburgh, Scotland, 2005.
- Morabito N, Crisafulli A, Vergara C, Gaudio A, Lasco A, Frisina N. Effects of genistein and hormone-replacement therapy on bone loss in early postmenopausal women: a randomized double-blind placebo-controlled study. *J Bone Miner Res*, 2002; 17(10):1904–12.
- Mulsow K, Eidenschink J, Melzig MF. FT-IR Method for the quantification of isoflavonol glycosides in nutritional supplements of soy (*Glycine max* (L.) Merr.). *Sci Pharm*, 2015; 83(2):377–86.
- Prabaningdyah NK, Riyanto S, Rohman A. Application of FTIR spectroscopy and multivariate calibration for analysis of curcuminoid in syrup formulation. *J Appl Pharm Sci*, 2018; 8(03):172–9.
- Rohman A. Analysis of curcuminoids in food and pharmaceutical products. *Int Food Res J*, 2012; 19:19–27.
- Rohman A, Ramadhani D, Nugroho A. Analysis of curcumin in *Curcuma longa* and *Curcuma xanthorrhiza* using FTIR spectroscopy and chemometrics. *Res J Med Plant*, 2015; 9:179–86; doi:10.3923/rjmp.2015.179.186
- Rohman A, Setyaningrum DL, Riyanto S. FTIR spectroscopy combined with partial least square for analysis of red fruit oil in ternary mixture system. *Int J Spectroscopy*, 2014; 1-5.
- Shao S, Duncan AM, Yang R, Marcone MF, Rajcan I, Tsao R. Systematic evaluation of pre-HPLC sample processing methods on total and individual isoflavones in soybeans and soy products. *Food Res Int*, 2011;44(8) 2425-2434.
- Sim C, Hamdan M, Ismail Z. Assessment of herbal medicines by chemometrics-assisted interpretation of FTIR spectra. *Anal Chim Acta*, 2004; 570: 1–14.
- Sulistiyowati E, Martono S, Riyanto S, Lukitaningsih E. Development and validation for free aglycones daidzein and genistein in soybeans (*Glycine max* (L.) Merr.) using RP HPLC method. *Int J Appl Pharm*, 2019; 11(2): 138-142.
- Wulandari R, Sudjadi, Martono S, Rohman A. Liquid chromatography and Fourier transform infrared spectroscopy for quantitative analysis of individual and total curcuminoid in *Curcuma longa* extract. *J App Pharm Sci*, 2018; 8(09):107–13.
- Yatsu FKJ, Koester LS, Bassani VL. Isoflavone-aglycone fraction from *Glycine max*: a promising raw material for isoflavone-based pharmaceutical or nutraceutical products. *Rev Bras Farmacogn*, 2016; 26(2):259–67.

**How to cite this article:**

Sulistiyowati E, Martono S, Riyanto S, Lukitaningsih E, Rohman A. Rapid quantitative analysis of daidzein and genistein in soybeans (*Glycine max* (L.) Merr.) using FTIR spectroscopy and multivariate calibration. *J Appl Pharm Sci*, 2020; 10(11):117–123.